



High Sensitivity Tools for Vaccine Development and Quality

Kevin Dahl, PhD and Linda Kidder, PhD

Hosted by Julie Nguyen and HORIBA Scientific



particlesell.com

Kevin Dahl, PhD - 20 years of particle and spectroscopic experience in pharma



Technology Consulting

- Instrumentation, laboratory



Data Consulting

- Reanalysis, interpretation, CMC support



Document Consulting

- Drafting and Review



Method Consulting

- Development, Optimization, Qualification, Validation



Experimental Consulting and Services

- Study design, execution, analysis, interpretation



Training Services

- Selected Instruments



kd@particlesell.com



- ◆ *The reward for work well done is the opportunity to do more.*
- Jonas Salk



Challenges in Vaccine Development

- ◆ Discovery
 - ◇ Rapid viral evolution – antigenic shift
- ◆ Formulation
- ◆ Storage and transport
 - ◇ Freeze-thaw
- ◆ Stability (ICH Q5C)
 - ◇ Toxic/immunogenic species
- ◆ In-use stability
 - ◇ Thaw/recon -> patient
- ◆ Release testing
 - ◇ Quality (ICH Q6B/Q6A)

**Particulate and
Spectroscopic
Analytical Tools**





Regulatory

- Modern vaccines inhabit a gray area between biotherapeutics (DS) and small molecule suspensions (DP)
 - Drug substance covered by ICH Q5A(R1) (safety)
 - Above deal with viral contamination/clearance, concentration, etc.
- What about drug product characteristics as CQAs?
- WHO Expert Committee on Biological Standardization, DRAFT Evaluation of the quality, safety and efficacy of messenger RNA vaccines for the prevention of infectious diseases: regulatory considerations, 2021

Particle size distribution (purity, consistency, safety)

light scattering such as dynamic or static light scattering; nanoparticle tracking analysis; electron microscopy; size-exclusion chromatography



Vaccines and Analytics


● Vaccine Delivery

- Live Attenuated or Inactivated Virus - Polio
- Adjuvanted - Nuvaxovid
- Lipid Nanoparticle (LNP) – Pfizer/Moderna
- Viral Vector/AAV – J&J
- Virus Like Particle (VLP) – Gardasil

} **PARTICLES!**


● Relevant Size Ranges and Common Techniques

- Cryo-TEM
- DLS
- NTA
- Laser Diffraction (LD)





Nanoparticle Tracking



- Based on video taken of point-scattering from particles

- Particle movement in suspension described by Brownian motion:

$$D_h = \frac{k_B T}{3\pi\eta D_m}$$

- Light scattering techniques measure hydrodynamic radius (D_h)

- Strengths


- Polydisperse (aggregated) samples
- Single-particle (number-based) technique
- Less sensitive to oversized material (relative to DLS)
- Governed by ISO 19430, ASTM E2834-12



Experimental: Vaccine Analytics


- Lyophilized sample with included diluent (1 mL)
 - Modified live virus with adjuvant (proprietary), sub-Q delivery, pH 6-8
 - Reconstitution after aseptic diluent addition, “SHAKE WELL”
- Sample diluted 10^5 in PBS
 - Concentration range match for NTA instrument
- 3-laser NTA system (445/520/635 nm)
 - 350 μ L sample stirred in cuvette
 - Practical limit of ~ 50 nm
 - Fluorescence capable

Reference: <https://www.durvet.com/product/canine-spectra-6/>






How things started:



- DLS results show a broad peak centered at 200 nm
 - Width suggests polydisperse sample – not optimal for DLS
 - Resolving power is $\sim x3/x5$ for two modes
 - Intensity-based distribution – r^6 intensity dependence
- Number-based NTA shows a single mode at 100 nm
 - DLS signal dominated by a small amount of larger debris
 - Reconstitution/dilution?
 - NTA signal characterizes main population in sample



NTA Results: Dilution Only




- 40 min collection time for six replicates
 - 50% decrease in particle concentration after 40 min
 - Access to concentration is critical!
 - D90 of 140 nm (no prominent aggregation observed)
- Change in sample indicates instability
 - Dilution in PBS destabilizes the formulation

Sample	Rep	Counts	D50 (nm)	D90 (nm)	Conc. (part/mL)
Spectra 6 - PBS	1	5349	98.9	136.3	1.0E+13
Spectra 6 - PBS	2	4531	99.9	139.2	8.6E+12
Spectra 6 - PBS	3	4067	99.9	140.0	7.8E+12
Spectra 6 - PBS	4	3355	100.6	137.9	6.4E+12
Spectra 6 - PBS	5	3134	99.9	139.5	6.0E+12
Spectra 6 - PBS	6	2730	100.0	142.9	5.2E+12
Average		3861	99.9	139.3	7.4E+12
%RSD		25.3	0.5	1.6	25.3



NTA Results: Dilution w/ 0.02% PS20




- Addition of 0.02% PS20 with PBS
 - PS20 should provide surface protection of unstable charge
- 40 min collection time for six replicates
 - ~15% decrease in particle concentration after 40 min
 - D90 has dropped to 132 nm

Sample	Rep	Counts	D50 (nm)	D90 (nm)	Conc. (part/mL)
Spectra 6 - PBS+0.02% PS20	1	5981	97.9	134.4	2.3E+13
Spectra 6 - PBS+0.02% PS20	2	5786	97.4	133.6	2.2E+13
Spectra 6 - PBS+0.02% PS20	3	5762	96.7	133.0	2.2E+13
Spectra 6 - PBS+0.02% PS20	4	5123	96.3	131.7	2.0E+13
Spectra 6 - PBS+0.02% PS20	5	5121	94.8	129.8	2.0E+13
Spectra 6 - PBS+0.02% PS20	6	5003	94.5	131.2	1.9E+13
	Average	5463	96.3	132.3	2.1E+13
	%RSD	7.79	1.44	1.29	7.79



NTA Results: Dilution w/ 0.2% PS20




- Addition of 0.2% PS20 with PBS
- 40 min collection time for six replicates
 - No decrease in particle concentration after 40 min
 - D90 135 nm
 - PS20 reversing “hidden” aggregation?
 - Distribution width shifting and narrowing

Sample	Rep	Counts	D50 (nm)	D90 (nm)	Conc. (part/mL)
Spectra 6 - PBS+0.2% PS20	1	7665	104.1	143.9	2.9E+13
Spectra 6 - PBS+0.2% PS20	2	7852	100.6	140.6	3.0E+13
Spectra 6 - PBS+0.2% PS20	3	8009	97.2	135.2	3.1E+13
Spectra 6 - PBS+0.2% PS20	4	7599	96.0	133.5	2.9E+13
Spectra 6 - PBS+0.2% PS20	5	7904	93.8	129.8	3.0E+13
Spectra 6 - PBS+0.2% PS20	6	7100	94.4	129.6	2.7E+13
	Average	7688	97.7	135.4	2.9E+13
	%RSD	4.24	4.07	4.28	4.24



NTA Results: Additional Work




- Three stress conditions tested for 48 hours
 - Shaker plate, direct sunlight, 50 °C
- Sample tolerates shaking
- Sample tolerates elevated temperature well
- Sample does not tolerate direct sunlight (UV)

Sample	Cond.	Counts	D50 (nm)	D90 (nm)	Conc. (part/mL)
Spectra 6 - PBS	None	3861	99.9	139.3	7.4E+12
Spectra 6 - PBS+0.2% PS20	None	7688	97.7	135.4	2.9E+13
Spectra 6 - PBS+0.2% PS20	48 hr. Shake	6693	98.9	133.6	2.6E+13
Spectra 6 - PBS+0.2% PS20	48 hr. Sun	3648	93.6	139.6	1.4E+13
Spectra 6 - PBS+0.2% PS20	50C	6828	91.2	137.4	2.6E+13



In Addition...



- Laser diffraction records scattering angle vs. intensity of scatter from suspension
 - Widely used for powder/high conc. suspensions
 - Minimum volume ~5 mL
 - Volume-weighted (r^3) distribution, can go below 100 nm
- LD is applicable to PS determination in vaccines
 - Products contain high concentration of particles
- Spectra 6 product shows band at 10 μm in addition to adjuvant band
 - Likely due to incomplete reconstitution
 - Explains poor DLS performance for same product!



Summary

- Vaccines are analytically 'different'
 - ◊ Biologic product necessitates specific testing
 - ◊ Particles are the product, testing with new technologies
- Incorporation of particle techniques
 - ◊ Directly measure CQA of Drug Product
 - ◊ Useful in formdev, stability, transport, and in-use studies
 - ◊ Can be readily worked into product release workflow
- Kevin Dahl, PhD: kd@particlesellc.com